

evolution^{neo+}

Clinical Chemistry Analyzer



Let our experience....
be your guide....

CONTACT US

+91 94 1683 1111



E-MAIL US

info@serachem-diagnostic.co.in



VISIT US

www.serachem-diagnostic.co.in

Technical Parameters

evolution^{neo+}

7 inches Color Touch screen interface
Large Graphic display with real time graph display for enzymatic reaction
Graph zoom in option providing clarity of the reaction curves
Built-in thermal printer to print results, enzymatic curves and standard calibration
curves India's First cloud based biochemistry analyzer,
Reagents management alarm, On board Incubation.

Multi-Programming Modes Absorbance
Endpoint with standard, factor, sample blank and
differential Fixed Time(Two-point) and multi point Kinetics
Multi standard calibration
Mono or bi-chromatic reading Turbidimetry

Optical System
Principle: photometry
Range of measurement: 0-2.5 A for 340nm;
0-3.0 A for 405nm onwards
Spectral range selection: automatic via filter wheel
Filters: 340, 405, 505, 546, 578 and 630nm plus two free
position Bandwidth: 10nm (+/- 2nm)
Light source: Halogen lamp 12 V20W
Blanking: Automatic
Detector: Silicon Photodiode

Memory
150+ User Defined Test
Programs 2000 Test Result
Reagent Blank
3 Level of QC

Power/Dimensions
Voltage: 120-230 VAC
Frequency: 50Hz
Maximum power: < 100 VA
Maximum relative humidity:
75% Compact form factor
Weight: ~5kgs

User Interface
Graphic LCD back Lit. 240x64 pixels
Touch Screen input
Data port
Type: Graphic thermal printer
Characters perline: 384 dots/line
Paper: Thermosensitive - 57mm
wide

Multi-Programming Modes Absorbance
Endpoint with standard, factor, sample blank and
differential Fixed Time(Two-point) and multi point Kinetics
Multi standard calibration
Mono or bi-chromatic reading Turbidimetry

Photometric Block
Peltier regulated 25 C, 30 C and 37 C
Accuracy: +/- 0.1 C
Flow cell of 33 u L
Removable cuvettes: semi-micro and Micro
Continuous flow system with perstatic pump
incorporated Programmable Sipping Volume

Quality Control Chart
Levy-Jennings chart using westgard's rules
30 Normal and abnormal QC values for selected
test



* with dedicated system packs

“Let our experience...
be your guide...”

www.serachem-diagnostic.co.in



EC REP

Serachem Diagnostic Limited Serachem Diagnostic, LLC
71-75, Shelton Street, Convent # 354 E, Main Street Suite, 1011 Middletown
Garden, LONDON, WC2H9JQ, D.E. 19709, USA, PH: +1 (302) 601-7242
ENGLAND

E-mail: serachem@serachem-diagnostic.com
www.serachem-diagnostic.com

Manufactured By

Serachem Diagnostic India Private Limited
82, Ward No. 10, Vashisth Colony, Near OHM
Convent School, Kurukshetra, Haryana-136118
E-mail: info@serachem-diagnostic.co.in
PH : +91-9416831111
www.serachem-diagnostic.co.in

Manufacturing Licence No.: MFG/IVD/2023/000121

Copyright © Serachem Diagnostic India Private Limited